



LATEX SURGICAL GLOVES

Description:

Sizes: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
Type: Powdered or Powder-free
Material: 100% Natural high quality latex
Color: Milky White

Key Feature: Soft; Comfortable naturally fits;
Textured Surface; Beaded cuff; Low level proteins;
Leak resistance; Tear resistance
Sterilization: Sterilized by Ethylene Oxide/ Gamma ray

Application: Healthcare; Surgical Room and Medical Room; Single use only
To help prevent contamination between caregivers and patients

Package: 1pair/bag, 50pairs/box, 10boxes/carton
Carton size: 67*26*40cm



LATEX EXAMINATION GLOVES

Description:

Size: S M L XL
Type: Powdered or Powder-free
Material: 100% Pure Natural latex
Color: Milky White

Key Feature: Soft; Comfortable naturally fits;
Textured Surface; Beaded cuff; Low level proteins;
Leak resistance; Tear resistance

Sterilization: Sterilized by Ethylene Oxide/ Gamma ray
Application: Healthcare; Homecare; Food Service; Laboratory/Industrial Safety;
To help prevent contamination between caregivers and patients

Package: 100pieces/box, 10boxes/carton
Carton size: 34*25.8*23cm



VINYL EXAMINATION GLOVES

Description:

Size: S M L XL
Type: Powdered or Powder-free
Material: High Quality PVC Material
Color: Clear, Red, Yellow, Green, Black, Blue; Customized

Key Feature: Soft; Comfortable naturally fits;
Compound material; Environmental; Good feeling;
Solid and durable
Application: Healthcare; Homecare; Food Service; Laboratory/Industrial Safety;
To help prevent contamination between caregivers and patients

Package: 100pieces/box, 10boxes/carton
Carton size: 34*25.8*23cm



FACE MASK EAR LOOP


Description:

Face Masks 3 Ply with Ear loops

- * This classical style Face Masks 3 Ply with Ear loops
- * with easy elastic ear loop and nose bar, 3 layer pleated design
- * Color: White/Blue/Green/Yellow/Pink/White
- * Breathable, soft lining, latex and fiberglass free

Non woven, 3Ply, with or without melt-blown fabric, ear loop. 17.5cm*9.5cm, 20+20+25g, light blue or white color, 50pc/box, 40boxes/carton.

The face masks are made of special non woven materials. It takes into consideration both the comfort of human body and the basic requirements of efficient breath. A flexible nose bar is adopted at the bridge of nose; the four sides of the mask join flatly and tightly to the face; hence, it avoids unflinching air directly entering into the mask and thus make breathe safely. Furthermore, her ear loop which is made of antibacterial spandex prevents any flu danger.



FACE MASK WITH TIE

Description:

Non woven, 3Ply, with or without melt-blown fabric, face mask with tie. 17.5cm*9.5cm, 20+20+25g, light blue or white color, 50pc/box, 40boxes/carton.

The face masks are made of special non woven materials. It takes into consideration both the comfort of human body and the basic requirements of efficient breath. A flexible nose bar is adopted at the bridge of nose; the four sides of the mask join flatly and tightly to the face; hence, it avoids unflinching air directly entering into the mask and thus make breathe safely. Furthermore, the ear loop which is made of antibacterial spandex prevents any flu danger.



SURGICAL CAP

Description:

Size: 61x13cm, 62x14cm; Customized
Weight: 20~40g/ m²
Type: With Ties
Material: Polypropylene, SMS Material
Color: Blue, Green, White; Customized

Key Feature: Comfortable naturally;
Compound material; Environmental;
Breathable; Solid and durable

Application: Healthcare; Homecare;
Food Service; Laboratory/Industrial Service;
To prevent dust and hairs falling

Package: 100pieces/pack, 10packs/ctn
Carton size: 28*28*30cm



BOUFFANT CAP (ROUNDED)

Description:

Elastic edge secure for fit
Light, soft, comfortable to wear
With effective protection
Economical, elegant
Regular: PP 10g/12g/18g/23g/28g/30g
Size: 18", 20", 21", 24" or customized
Color: white, blue, green
Package: 100pcs/bag, 1000pcs/ctn



BOUFFANT CAP (CLIP CAP/MOP CAP)

Description:

Light, soft, high air permeability
Machine manufactured, no manual handling
Economical yet effective headwear
Single or double elastic edge can be selected
Comfortable for long durations
Regular: PP 10g /12 g
Size: 19" 20" 21"
Color: white, blue, green, red
Package: 100pcs/bag, 2000 pcs/ctn



SURGICAL GOWN

Description:

Size: 115×127cm, 115×136cm,
115×145cm, 130×150cm; Customized
Gram Weight: 40gsm (30-55gsm are available)
Material: PP, PE+PP, SMS or Spunlace
non-woven fabric
Color: Blue; Customized

Key Feature: Soft and light weight material
provides long term comfort for long procedures;
Water-repellent, oil-repellent, alcohol repellent,
blood repellent and antistatic

Application: Healthcare; Surgical Room and
Medical Room; Chemical industrial use
To help prevent contamination between
caregivers and patients

Package: 2000pieces/carton
Carton size: 40*25*35cm



SURGICAL SHOES COVER


Description:

Size: 15x41cm/15x36cm/15x39cm; Customized
Weight: 1.2g~4.0g/piece
Material: Non-woven Material
Color: Blue, Green, Clear, Red, Yellow; Customized

Key Feature: Comfortable naturally; Compound material;
Environmental; Breathable; Solid and durable

Application: Healthcare; Homecare; Food Service;
Laboratory/Industrial Service;
To help filter pollen, dust and bacteria

Package: 2000pieces/carton
Carton size: 42*28*52cm



FIRST AID BAND

Description:

First aid band, first aid, first aid products
PE, PVC, PU, POE, elastic cotton, cotton,
non woven

Size: 70mmx18mm, 72mmx19mm,
75mmx22mm, Ø22mm, cross type,
butterfly type, H type etc



ALCOHOL PAD

Description:

3*6cm pad with 5*5cm pouch, 70% isopropyl alcohol, individually sealed swabs, 100pcs/box, 200 pcs/box, 10000pcs/carton
Alcohol swab, with 70% IPA.

Original solution from Exxon Mobil 99.8% solution.
60x30mm, single pack.

Standard packing:

English neutral packing, 100pc/box or 200pc/box, 100box/carton;
or as client requirement.

Specs:

Pouch size: 50*50±2mm

Pouch Material: Aluminum foil paper, 73g/m2

Non-woven size: 60*30±2mm,

Solution: 70% isopropyl alcohol

Solution Weight: 0.4g- 0.5g

The material of the box: cardboard with white face and mottled back; 300g/m2;



POVIDONE-IODINE PREP PAD

Description:

One 3*6cm Prep Pad in 5*5cm pouch saturated with a 10% Povidone Iodine Solution equivalent to 1% available iodine. 100pc/box, 10000pc/carton. English neutral packing

Saturated with 10% Povidone- Iodine, solution equivalent to 1% Povidone- Iodine

Pouch Material: Aluminum foil paper, 90g/m2

Non-woven size: 60*30±2mm,

Solution: with 10% Povidone- Iodine, solution equivalent to 1% Povidone- Iodine

Solution Weight: 0.4g- 0.5g

The material of the box: cardboard with white face and mottled back; 300g/m2;



SILK SURGICAL TAPE

white colour, skin colour

Description:

Length: 5yard, 10yard
Width: 1.25, 2.5, 5, 7.5, 10cm, etc

Features:

1. With reliable adhesiveness, low sensitization, excellent compliance, no residue glue, no harm to skin and air permeability, it enables the skin to breath freely;
2. Most are easy-to-tear products, making usage very convenient and comfortable;

Application

- *Apply to clean and dry skin;
- * When adhesive, please gently press the adhesive tape from center to the two sides, Making it closely combine the skin;
- * Please don't adhere the center of adhesive tape to the wound.



ZINC OXIDE PLASTER


Description:

Length: 5yard, 10yard
Width: 1.25, 2.5, 5, 7.5, 10cm, etc

- 1) 100%cotton adhesive bandage
- 2) transparent plastic. Adhesive bandage
- 3) non-wound adhesive plaster

Feature:

- 1) Good adhesive ability
- 2) Easy to store, long storage life
- 3) Harmless to skin



NON-WOVEN ADHESIVE TAPE

white colour, skin colour

Description:

Length: 5yard, 10yard
Width: 1.25, 2.5, 5, 7.5, 10cm, etc

It can fix and keep wound dressings, bandages, catheters and replacing medical materials. Permeability makes the skin breathe naturally, and excrete water, air and sweat, and also reduces the possibility of infection. Adaptation gentle and thin fiber can make the patients feel comfortable and suit the curve of people's body. Hypoallergenic and proper viscosity will not irritate the skin and be firmly adhered. Do not hurt the skin when it released from the skin.



MICROPORE PE ADHESIVE TAPE

Transparent

Description:

Length: 5yard, 10yard
Width: 1.25, 2.5, 5, 7.5, 10cm, etc

The substrate is perforated PE, moisture and sweat can be drained out easily. The initial sticking force of this adhesive tape is enough, no residue on the skin.

Features:

1. PE material with micro-holes
2. Hypoallergenic
3. Good air permeability
4. Soft and comfortable
5. Easy tear by hand in longitude and latitude directions



COTTON WOOL

Description:

50g, 100g,200g, 250g, 400g, 454g, 500g, 1000g

Medical defatted absorbent cotton fibers from 100% cotton, produced two oxygen bleach, the whiteness of 85 degrees. With strong absorption capacity for beauty, cosmetics, cleaning trauma wounds and absorb the blood alcohol disinfection cleaning of medical equipment and instruments.



PLASTER OF PARIS BANDAGE

Description:

Size: 5cm x 2.7m, 7.5cm x 2.7m, 10cm x 2.7m, 15cm x 2.7m, 20cm x 2.7m, 25cm x 2.7m, 30cm x 2.7m. Etc.
Various size is available.

Consists of plain-weave gauze fabric coated with blend of alpha and beta calcium caliphate crystals, spooled on round perforated plastic core.

Benefits

- Immersion time 2 to 3 seconds only.
- Excellent moulding ability.
- Initial setting time within 3 to 5 minutes, at immersion water temperature of 20 °C.
- Can be carefully subjected to weight bearing after 30 minutes.
- Very low plaster loss.

- When completely hardened have high strength at low bandage consumption.

Indications

- Immobilization of all types of fractures.
- Immobilization of diseased bones and joints.
- Correction of deformities.



NON WOVEN COHESIVE FLEXIBLE BANDAGE

It is an economical, hand tear, self-adhering elastic bandage. It is soft, easy applied, sticks only to itself and allows skin to breath. It can be used in many of the situations where you could use an elastic bandage. It is lighter and more conforming to your shape.

Features:

- 1.Elastic non woven substrate
- 2.Sticks to itself, does not stick to hair, skin, clothing, no pins and clips needed
- 3.Soft, breathable and comfortable
- 4.Easy hand tear, no scissors needed
- 5.Provide light compression, apply properly to avoid cutting circulation
- 6.Stable and reliable cohesiveness
- 7.Constant unwinding tension
- 8.Good tensile strength

Size	Color	Packing(rolls/box)
2.5cm x 450cm	White/Tan/ Red /Dark blue/ Green/ Yellow/ Neon color	72
5.0cm x 450cm		36
7.5cm x 450cm		24
10cm x 450cm		18
15cm x 450cm		12



LATEX FREE NON WOVEN COHESIVE FLEXIBLE BANDAGE

It is contains no natural rubber latex.



Features:

- 1.Elastic non woven substrate
- 2.Sticks to itself, does not stick to hair, skin, clothing, no pins and clips needed
- 3.Latex-free, causes no latex induced allergic reactions
- 4.Soft, breathable and comfortable
- 5.Easy hand tear, no scissors needed
- 6.Provide light compression, apply properly to avoid cutting circulation
- 7.Stable and reliable cohesiveness
- 8.Good tensile strength
- 9.Water resistant

Size	Color	Packing(rolls/box)
2.5cm x 450cm	White/Tan/ Red /Dark blue/ Green/ Yellow/ Neon color	72
5.0cm x 450cm		36
7.5cm x 450cm		24
10cm x 450cm		18
15cm x 450cm		12



COTTON COHESIVE FLEXIBLE BANDAGE

It is a kind of cotton cohesive flexible bandage. Its cotton substrate is soft, breathable and easy to tear. The cohesive property allows it to stick to itself but not to skin or hair.

Features:

- 1.Elastic cotton substrate
- 2.Sticks to itself, does not stick to hair, skin, clothing, no pins and clips needed
- 3.Soft, breathable and comfortable
- 4.Easy hand tear, no scissors needed
- 5.Provide light compression, apply properly to avoid Cutting circulation
- 6.Stable and reliable cohesiveness
- 7.Constant unwinding tension
- 8.Good tensile strength
- 9.Competitively priced

Size	Color	Packing(rolls/box)
2.5cm x 450cm	White/Tan/Black	72
5.0cm x 450cm		36
7.5cm x 450cm		24
10cm x 450cm		18
15cm x 450cm		12



LATEX FREE COTTON COHESIVE ELASTIC BANDAGE

It is a latex free cotton cohesive flexible bandage. It sticks to itself but not to hair or skin. It does not cause latex induced allergic reactions.



Features:

- 1.Elastic cotton substrate
- 2.Sticks to itself, does not stick to hair, skin, clothing, no pins and clips needed
- 3.Latex-free, causes no latex induced allergic reactions
- 4.Soft, breathable and comfortable
- 5.Easy hand tear, no scissors needed
- 6.Provide light compression, apply properly to avoid cutting circulation
- 7.Stable and reliable cohesiveness
- 8.Good tensile strength

Size	Color	Packing(rolls/box)
2.5cm x 450cm	White/Tan/Black	72
5.0cm x 450cm		36
7.5cm x 450cm		24
10cm x 450cm		18
15cm x 450cm		12



COTTON ELASTIC ADHESIVE BANDAGE

It is a heavy duty elastic adhesive bandage with aggressive adhesion. Its soft cotton have very good conformability. It is applied when moderate to maximum strength is necessary.



Features:

- 1.Cotton elastic substrate
- 2.Latex-free, causes no latex induced allergic reactions
- 3.Soft and comfortable
- 4.Heavy and stable stretch
- 5.Provide moderate to maximum compression, apply properly to avoid cutting circulation
- 6.Strong and reliable stickiness
- 7.Constant unwinding tension
- 8.Leaves no residue on body parts
- 9.Color thread in the middle of the bandage facilitating overlapping
- 10.Reverse coating is available

Size	Color	Packing(rolls/box)
2.5cm x 450cm	White/Skin	48
5.0cm x 450cm	White with yellow middle line/ Skin with red middle line	24
7.5cm x 450cm		16
10cm x 450cm		12



HIGH ELASTIC BANDAGE

It can be wrapped and secured easily. It provides patient comfort and is durable enough to be reusable

Features:

- 1.Latex-Free or not
- 2.Double closure delivers maximum control and consistent compression
- 3.Unique knitted construction
- 4.Superior absorbency and breath ability soft and comfortable
- 5.Require no pins, clip or tapes

Size	Color	Packing(rolls/box)
5.0cm x 450cm	Natural	10
7.5cm x 450cm		10
10cm x 450cm		10
15cm x 450cm		10



CREPE BANDAGE

Features:

- 1.The product has the features of good stretching, strong permeability, uniform pressure, snug and fitting binding, no slippage etc.
- 2.It has the auxiliary treatment function to strain of limbs, injury of soft tissue and arthrodynia etc, and can also the used as splints, varicose vein under supplement treatment and limbs fracture etc.

Size	Color
2.5cm x 450cm	White/Natural
5.0cm x 450cm	
7.5cm x 450cm	
10cm x 450cm	
15cm x 450cm	



ORTHOPEDIC CASTING BANDAGE



This product is made from high active polyurethane resin and fabric, it has good biocompatibility; the experiments of animals acute and chronic toxicity test confirmed that this product is non-toxic and non-distortion. It mainly used in orthopedics, plastic surgery and general surgery for fractures and sprains of the external fixation and to promote fracture healing. This is the new substitute of the traditional plaster bandage.

Features:

1. High hardness, light weight
2. Curing speed
3. Excellent X-ray projection
4. Good waterproof
5. Clean, Sanitary

Size	Color	Packing
5.0cm x 360cm	White/Green/ Blue/ Pink/Yellow/ Orange	10roll x 12box
7.5cm x 360cm		10roll x 12box
10.0cm x 360cm		10roll x 9box
12.5cm x 360cm		10roll x 9box
15.0cm x 360cm		10roll x 9box

SURGICAL FILM



Features

- 1)Transparency: The film is ultra thin and high transparent, hence clear observation to surgical operation part.
- 2)Permissibility: No effect to normal breath of the skin, protection from water buildingup under the film, thus reducing infection and providing a germ-free environment for surgical cut.
- 3)Water & Germ Resistance: Prevention from invasion of water, germ or foul and restrain wandering about of germ during operation.
- 4)High elasticity: Compliant binding on skin along body curve, hence reliable binding.
- 5)Low sensitivity & moderate viscosity: no irritation to wound and no injury to skin when removing the dressing.

Description for use

Protection of operation parts before and during the operation and prevention from cross infection.

Direction for use


- 1)Disinfect the operation parts, then thoroughly dry the operation field area.
- 2)Spread the film to extend by two sides, peel the liner and put the film on operation area.
- 3)Smooth the film from the center to circumference; peel the PE film, then operation.
- 4)After operation lifting the film near the wound and sew up. The last, remove the surgical film.
- 5)Pay attention not to remove together the conduit or other devices covered by it.

Specifications

Code	Specification	Packing
AP1420	14cm X 20cm	20
AP2030	20cm X 30cm	20
AP3030	30cm X 30cm	20
AP3045	30cm X 45cm	20
AP4545	45cm X 45cm	20
AP4560	45cm X 60cm	10
AP4565	45cm X 65cm	10
APE2030	20cm X 30cm	20
APE3030	30cm X 30cm	20
APE3045	30cm X 45cm	20
APE4545	45cm X 45cm	20
APE4560	45cm X 60cm	10
APE4565	45cm X 65cm	10

The products can be made on customer' request.

ADHESIVE WOUND DRESSING



Features

- 1)Permlssliblity: the network structure of spun-laced nonwoven cloth can make skin breathe freely, get rid of steam and sweat, thus reducing the occurrence of wound infection effectively.
- 2)Low sensitization: the medical pressure-sensitive adhesive to bind on skin has moderate viscosity: has no irritation to wound and no injury to skin when removing the dressing.
- 3)No binding to wound: the absorptive pad with network cover will not stick, and can absorb effusion effectively without any sticking to wound; thus causing no pain when removing the pad.
- 4)Compliance: being soft, light and elastic, the material can comply with body outline and curves without any obstacle to muscle action.

Description for use

Protects the wound site after surgical operations.

Direction for use

- 1)Prepare the wound according to institution protocol.Allow all cleansing solutions and skin peotectants to dry completely.
- 2)Peel the liner from the dressing, bind the dressing on wound and press the circumference to make it firm.


Caution

- 1)Forbidden for wound infected and ulcerate d. Please stop using it or adopt relevant means according to doctor's diagnosis, if any hyperemia, swelling, seepage, fever, et c. after infection.
- 2)It couldn't be used as fixing aim of artery conduit.
- 3)It can't replace sewing, anastasis, skin disinfection, drying, etc measures.
- 4)Please select the size & specification enough to assure its binding firmly on dry healthy skin around the wound when use, and no strain is made by stretching.
- 5)Pay attention not to remove together the conduit or other devices covered by it.


Specifications

Code	Specification	Packing
AP0607	6cm X 7cm	20
AP0910	9cm X 10cm	20
AP0915	9cm X 15cm	20
AP0920	9cm X 20cm	20
AP0925	9cm X 25cm	20
AP0930	9cm X 30cm	20
AP0935	9cm X 35cm	20
AP0940	9cm X 40cm	20
AP0945	9cm X 45cm	20
AP0950	9cm X 50cm	20
AP0612(EYE)	6cm X 12cm	20
AP1212	12cm X 12cm	20
AP1212-1	12cm X 12cm	20
AP1010	10cm X 10cm	20
AP1015	10cm X 15cm	20
AP1020	10cm X 20cm	20

Needle identification & Application



Needle Shape



PRODUCT CODE GUIDE

Code Guide


THREAD TYPE	CODE	THREAD LENGTH	CODE	THREAD DIAMETER	CODE	NEEDLE LENGTH	CODE	NEEDLE CURVATURE	CODE	CROSS SECTION	CODE
Polyglycolic Acid	PGA	15 cm	1	11/0	110	6.5	06	Straight	A	● round bodied	0
Rapide Polyglycolic acid	SPGA	20 cm	2	10/0	100	7	07				
Polyglactin 910	PGL	30 cm	3	9/0	90	8	08				
Rapide polyglactin 910	MPGL	45 cm	4	8/0	80	8.5	08				
Poliglecaprone 25	PGLP	60 cm	5	7/0	70	10	10				
Polydioxanone	PDS	70 cm	6	6/0	60	11	11	1/2 circle	B	● round bodied heavy	1
Chromic Catgut	CC	75 cm	7	5/0	50	12	12				
Plain Catgut	CP	90 cm	8	4/0	40	13	13	1/2 circle (double)	C	▲ (curved)cutting	2
Silk (Braided)	SK	100 cm	9	3/0	30	14	14				
Polyester(Braided)	PB	120 cm	10	2/0	20	15	15	1/4 circle	D	▲ (curved)cutting heavy	3
Nylon monofilament (Polyamide)	NL	125 cm	11	1/0	10	16	16				
Polypropylene monofilament	PPM	150 cm	12			17	17	1/4 circle (double)	E	▼ reverse cutting	4
Stainless Steel Wire	SW	180 cm	13			18	18				
						19	19	3/8 circle	F	● tapercut	6
						20	20				
						21	21	3/8 circle (double)	G	● tapercut heavy	7
						22	22				
						23	23	5/8 circle	H	▼ micro-point spatula curved	8
						24	24				
						25	25	5/8 circle (double)	I	○ blunt taper point	9
						26	26				
						27	27	LOOP	J	○ blunt taper point heavy	10
						28	28				
						29	29				
						30	30				
						31	31				
						32	32				
						33	33				
						34	34				
						35	35				
						36	36				
						37	37				
						38	38				
						39	39				
						40	40				
						41	41				
						42	42				
						43	43				
						44	44				
						45	45				
						46	46				
						47	47				
						48	48				
						49	49				
						50	50				
						51	51				
						52	52				
						53	53				
						54	54				
						55	55				
						56	56				
						57	57				
						58	58				
						59	59				
						60	60				
						61	61				
						62	62				
						63	63				
						64	64				
						65	65				
						66	66				
						67	67				
						68	68				
						69	69				
						70	70				
						71	71				
						72	72				
						73	73				
						74	74				
						75	75				
						76	76				
						77	77				
						78	78				
						79	79				
						80	80				
						81	81				
						82	82				
						83	83				
						84	84				
						85	85				
						86	86				
						87	87				
						88	88				
						89	89				
						90	90				

Example Code: SK 7 10 28 B 2

SK71028B2
Product code(Example)

SURGICAL SUTURES

Sutures Materials



Absorbable

- Polyglycolic acid
- Rapide Polyglycolic acid
- Polyglactin 910
- Rapide polyglactin 910
- Poliglecaprone 25
- Polydioxanone
- Chromic catgut
- Plain catgut

Non-absorbable

- Silk braided
- Polyester braided
- Nylon monofilament
- Polypropylene monofilament
- Stainless steel wire





DISPOSABLE SURGICAL BLADES

Description:

Disposable stainless steel surgical blades
Disposable carbon steel surgical blades
(Sterilization with Gamma Radiation)

We can process the products according to the samples and pack the products according to the fixed brands which are supplied by our customers.

Packing: 100 pcs/box, 5000pcs/ctn, 34*20*17cm



Disposable Surgical Scalpels



SURGICAL SCALPEL

Disposable Surgical Scalpels with plastic handle

Description:

Form of packing: packed by composition aluminium foil, aluminium foil/PVC plastics.

Way of sterilization: sterilized by Gamma Radiation.

Packing: 10pcs/box, 500pcs/ctn



DISPOSABLE SAFETY SURGICAL SCALPEL


Description:

The Index safety scalpel supports and enhances the safe passing procedures and protocols already practiced in operating rooms and theatres.


The scalpel features a strong temporary forward lock position to prevent uncontrolled retraction when contacting with bone during use, an additional notch to hold the blade more securely in the passing position and a permanent rear lock to be engaged prior to final disposal in a sharps container.

Packing: 10pcs/box, 5000pcs/ctn, 34*20*17cm


How to Use




Press button and do not push forward




Pull guard back to expose blade. Listen for audible click when fully open and ready for use.



Push guard forward to cover blade. Listen for audible click when fully closed.



Hold one or both sides of guard. Push forward and hear audible click when locked.



Dispose of the Index Safety Scalpel in a suitable container.



SPINAL NEEDLES COMMON POINT

Description:

Latex Free
Convenient blister packaging
Siliconized for smooth comfortable
Penetration
Transparent plastic hub

Package:50pcs/box,20box/carton,1000pcs/ctn,
51cm*42cm*25.5cm

Description

16G(1.60mm)X3 1/2"
17G(1.40mm)X3 1/2"
18G(1.20mm)X3 1/2"
19G(1.10mm)X3 1/2"
20G(0.90mm)X3 1/2"
21G(0.80mm)X3 1/2"
22G(0.70mm)X3 1/2"
23G(0.60mm)X3 1/2"
24G(0.55mm)X3 1/2"
25G(0.50mm)X3 1/2"
26G(0.45mm)X3 1/2"
27G(0.40mm)X3 1/2"

Hub Color

□ White
■ Purple
■ Pink
■ Cream
■ Yellow
■ DeepGreen
■ Black
■ DeepBlue
■ Violet
■ Orange
■ Brown
■ Grey



SPINAL NEEDLES PENCIL POINT

Description:

Latex Free
Convenient blister packaging
Siliconized for smooth comfortable
Penetration
Transparent plastic hub

Package:50pcs/box,20box/carton,1000pcs/ctn,
51cm*42cm*25.5cm

Description

16G(1.60mm)X3 1/2"
17G(1.40mm)X3 1/2"
18G(1.20mm)X3 1/2"
19G(1.10mm)X3 1/2"
20G(0.90mm)X3 1/2"
21G(0.80mm)X3 1/2"
22G(0.70mm)X3 1/2"
23G(0.60mm)X3 1/2"
24G(0.55mm)X3 1/2"
25G(0.50mm)X3 1/2"
26G(0.45mm)X3 1/2"
27G(0.40mm)X3 1/2"

Hub Color

□ White
■ Purple
■ Pink
■ Cream
■ Yellow
■ DeepGreen
■ Black
■ DeepBlue
■ Violet
■ Orange
■ Brown
■ Grey




EPIDURAL NEEDLES

Description:

Latex Free
Convenient blister packaging
Siliconized for smooth comfortable
Penetration
Transparent plastic hub
Used in epidural puncture for continuous epidural
anesthesia.

Package:50pcs/box,20box/carton,1000pcs/ctn,
51cm*42cm*25cm



SINGLE USE BONE MARROW ASPIRATION NEEDLE

Description:

The needle edge is firm and sharp.
The positioning body may adjust the length of the needle tube.
The needle tube is scaled for the purpose of puncture depth.
The optional sample tube reduces the waste of sampling.

Package: 200pcs/ctn, 47.5cm*39.5cm*41cm

Size:	Length(mm)
1.2(18G)	25(10-25),40(25-40),55(40-55)
1.6(16G)	25(10-25),40(25-40),55(40-55)
1.8(15G)	25(10-25),40(25-40),55(40-55)
2.1(14G)	25(10-25),40(25-40),55(40-55)



SINGLE USE BONE MARROW BIOSPY NEEDLE

Description:

Used in bone puncture to scrub bone matrix and suck bone marrow
Variety edge of the needle meets the different clinical demand.
The optional connector of convertible may get 10mm or 15mm length of tissue biopsy.

Package:400pcs/ctn, 57.5cm*44.5cm*49.5cm

Size:	Length(mm)
2.4(13G)	70
2.7(12G)	70
3.0(11G)	100



BONE MARROW KIT

Description:

According to clinical demand, the company may provide suite of instruments and accessories of aspiration type, biopsy type, for children and combined aspiration/biopsy type for get bone marrow tissue for biopsy.

Single-use Bone Marrow Tray

	Name of Option	Amount	BMT-A	BMT-B	BMT-C	BMT-A/C
Basic Option	Basic Option Bone Marrow Aspiration Needle	1	*	-	-	*
	Bone Marrow Biopsy Needle	1	-	*	-	*
	Bone Marrow Needle for children	1	-	-	*	-
	Sampling tube	1	*	-	*	*
	Connecting column	1	-	*	-	*
Option for choice	sterile syringe, sterile injection needle, rubber gloves for medical use, prep sponge swabs, surgical tweezer, absorbent gauze, Infusion plaster, surgical towel and surgical drape etc.					



COMBINED SPINAL EPIDURAL ANESTHESIA KIT

Contents	Specifications
Epidural needle	One 16G/18G
Penpoint spinal needle	One 25G/27G L (one hole)
Epidural catheter	One pull-resistance
Epidural catheter adapter	One (Tuohy Borst)
Epidural catheter filter	One (green) 5 v
Epidural needle filter	One(yellow)0.5 v
LOR Glass syringe	One 5cc
Syringes	Two (3cc and 20cc)
Negative pressure tube	One
Fenestrated cover	One 60cm*60cm with hole 10cm
Bottom cover	One 50cm*50cm
Tray cover	One 70cm*95cm
Epidural catheter thread assist guide	One
Adhesive tapes	One 5cm*40cm, two 1.25*10cm
Epidural catheter adhesive slice	One
Compresses	Three 8layers 8cm*10cm
Spinal needle introducer	One 16G 30mm
Sponge brush	Three
Surgical latex gloves	A pair #7.5
Cannula	Three (OD 0.5mm L20mm,OD 0.7mm L32mm, OD 1.6mm L38mm)
Products tray	One



EPIDURAL ANESTHESIA KIT

Contents	Specifications
Epidural needle	One 15/16G/18G
Epidural catheter	One pull-resistance
Epidural catheter adapter	One (Tuohy Borst)
Epidural catheter filter	One (green) 5 v
Epidural needle filter	One(yellow)0.5 v
LOR Glass syringe	One 5cc
Syringe	One 20cc
Negative pressure tube	One
Fenestrated cover	One 60cm*60cm with hole 10cm
Bottom cover	One 50cm*50cm
Tray cover	One 70cm*95cm
Epidural catheter thread assist guide	One
Adhesive tapes	One 5cm*40cm, two 1.25*10cm
Epidural catheter adhesive slice	One
Compresses	Three 8layers 8cm*10cm
Spinal needle introducer	One 16G 30mm
Sponge brush	Three
Surgical latex gloves	A pair #7.5
Cannula	Three (OD 0.5mm L20mm,OD 0.7mm L32mm, OD 1.6mm L38mm)
Products tray	One



SPINAL ANESTHESIA KIT

Contents	Specifications
Spinal needle	One 20G/22G/25G
Epidural catheter filter	One (green) 5 v
LOR Glass syringe	One 5cc
Syringes	5cc
Fenestrated cover	One 60cm*60cm with hole 10cm
Bottom cover	One 50cm*50cm
Tray cover	One 70cm*95cm
Adhesive tapes	One 5cm*40cm, two 1.25*10cm
Compresses	Three 8layers 8cm*10cm
Spinal needle introducer	One 16G 30mm
Sponge brush	Three
Surgical latex gloves	A pair #7.5
Cannula	Three (OD 0.5mm L20mm,OD 0.7mm L32mm, OD 1.6mm L38mm)
Products tray	One



DISPOSABLE SYRINGES

Description:

High Quality Stainless Steel Needle
Larger Inner Diameter High Flow
Sharpness Maximize Patient Comfort
Color Coded Hub by Size for Clear Recognition
Semi-transparent Needle-Hub for Clearness of Flashback
Individual and Bulk Packaging, Sterile & Non-sterile

Luer Lock

Size	Lock Ref	Package
1ml	SY-1-L	3600pcs/ctn,65*38*49.5cm
2ml	SY-2-L	2400pcs/ctn,60*46*42cm
3ml	SY-3-L	2400pcs/ctn,60*46*42cm
5ml	SY-5-L	2000pcs/ctn,80*42*38cm
10ml	SY-10-L	1200pcs/ctn,56*46*42cm
20ml	SY-20-L	800pcs/ctn,58*38*52cm
30ml	SY-30-L	400pcs/ctn,47*38*43cm
50ml	SY-50-L	200pcs/ctn,45*33*44cm
60ml	SY-60-L	200pcs/ctn,45*33*44cm

Luer Slip

Size	Lock Ref	Package
1ml	SY-1-L	3600pcs/ctn,65*38*49.5cm
2ml	SY-2-L	2400pcs/ctn,60*46*42cm
3ml	SY-3-L	2400pcs/ctn,60*46*42cm
5ml	SY-5-L	2000pcs/ctn,80*42*38cm
10ml	SY-10-L	1200pcs/ctn,56*46*42cm
20ml	SY-20-L	800pcs/ctn,58*38*52cm
30ml	SY-30-L	400pcs/ctn,47*38*43cm
50ml	SY-50-L	200pcs/ctn,45*33*44cm
60ml	SY-60-L	200pcs/ctn,45*33*44cm



SYRINGE WITH CATHETER TIP

Description:

PE Or Blister Package
Sterilized By EO Gas, Non-toxic, Non-pyrogenic,
Single Use Only

Size: 50ml/60ml

Package: 5000pcs/ctn,45*45*28.5cm



AUTO DISABLE SYRINGES

Description:

1. Sterilized by EO gas, non-toxic, non-pyrogenic, single use only.
2. Prevent the re-use of syringes and thus the transmission of diseases.
3. The integrated needle is an important feature of the non-dead-space design, which reduces waste of medication to a minimum.
4. Equipped with the Terumo lancet point needle which guarantees unsurpassed needle quality and minimal patient trauma.

Size	Slip/Lock	Package
0.5ml	SY-0.5-L	2400pcs/ctn,57*38*35cm
1ml	SY-1-L	2400pcs/ctn,57*38*35cm
2ml	SY-2-L	2400pcs/ctn,51*41.5*32cm
3ml	SY-3-L	2000pcs/ctn,51*41.5*32cm
5ml	SY-5-L	1200pcs/ctn,53*41.5*32cm
10ml	SY-10-L	1200pcs/ctn,65*41.5*42.5cm



SAFETY DISABLE SYRINGE

Description:

1. Sterilized by EO gas, non-toxic, non-pyrogenic, single use only.
2. Prevent the re-use of syringes and thus the transmission of diseases.
3. Needle will be retracted into the barrel to reach fully destruction.

Size	Lock Ref	Package
0.5ml	SY-2-L	2400pcs/ctn,57*38*35cm
1ml	SY-3-L	2400pcs/ctn,57*38*35cm
2ml	SY-5-L	2400pcs/ctn,51*42*32cm
3ml	SY-10-L	2000pcs/ctn,51*42*32cm
5ml	SY-20-L	1200pcs/ctn,53*42*32cm
10ml	SY-30-L	1200pcs/ctn,65*42*43cm
20ml	SY-50-L	300pcs/ctn,65*35*36cm

HYPODERMIC NEEDLE

Description:

High Quality Stainless Steel
Larger Inner Diameter High Flow
Sharpness Maximize Patient Comfort
Color Coded Hub by Size for Clear Recognition
Semi-transparent Needle-Hub for Clearness of Flashback
Individual and Bulk Packaging, Sterile & Non-sterile
Package: 100pcs/1 box, 2000pcs/carton, 56x25.5x37.5cm



Size	Hub Color
27Gx1/2"	Grey
26Gx1/2"	Beige
25Gx5/8"	Dark Blue
25Gx1"	Dark Blue
25Gx1 1/2"	Dark Blue
23Gx3/4"	Light Blue
23Gx1"	Light Blue
22Gx3/4"	Black
22Gx1"	Black
22Gx1 1/2"	Black
21Gx1"	Green
21Gx1 1/2"	Green
20Gx1"	Yellow
20Gx1 1/2"	Yellow
18Gx1"	Pink
18Gx1 1/2"	Pink
19Gx1"	Brown
27Gx1 1/2"	Grey
27Gx1 1/4"	Grey
20Gx3/4"	Yellow
19Gx1 1/2"	Brown
20Gx2"	Yellow



SCALP VEIN SET

Description:

Butterfly shaped wings for better handling and attachment to skin
Slip or luer set end matches well with different types of infusion sets
Siliconed cannula for smooth penetration
Tube length: 15-30cm
Gauge 18g - 26g available
ETO sterilized
Internationally accepted color-code
Manufacturing using automated polishing and lubrication processes
Strict quality control procedures ensure consistent performance

Package: 100pcs/1 box, 2000pcs/carton, 56x25.5x37.5cm

Pipe outside diameter

Gauge	Colour code	O.D.mm	Length mm
26G	Brown	0.45	5/8"
25G	Orange	0.5	3/4"
24G	Purple	0.55	3/4"
23G	Blue	0.6	1"
22G	Black	0.7	1"
21G	Dark Green	0.8	1 1/8"
20G	Yellow	0.9	1 1/8"
19G	Cream White	1.1	1 1/4"
18G	Pink	1.2	1 1/4"



0.3CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size:29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



1/2CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size:28GX1/2" / 29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



1CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size:28GX1/2" / 29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



INSULIN PEN NEEDLES

Description:

Optimum treatment with maximum simplicity
Sterile, Non-toxic, Non-pyrogenic
Double point system
Patient proof with triple beveling for easy and atraumatic penetration
Cartridge point sharpened for a perfect punch with no fragmentation
Compatible with all available pen models
Easy flow of medication without any impurities

Size:29G(0.33)X1/2"(12mm) / 31G(0.25)X1/4"(6mm)
/ 31G(0.25)X1/3"(8mm)

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



STERILE INFUSION SETS

Description:

COMMON FEATURES

Soft and kink resistant PVC tubing
Clear, transparent and flexible drip chamber
Moulded latex flush ball injection membrane for intermittent medication
Accurate flow controller
Approximately 15 drops / ml , 20 drops / ml, 60 drops / ml

Tube Length : 150 cm , 160 cm , 180 cm , 200 cm
Tube : I.D. : 3.0 mm , O.D. : 4.1 mm

Packing:PE package or paper-plastic package

Options available:

- With or without air vented spike
- With or without needle
- With or without "Y" injection port (Latex or Latex - free)
- Injection site material with latex or polyisoprene
- Luer lock or luer slip connector
- With or without flow regulator integrated into Infusion set for precise flow control
- With or without filter

Size:18G/19G/20G/21G/22G/23G/24G/25G/26G/27G

Package:500pcs/ctn, 58*38*36cm



BLOOD TRANSFUSION SETS

Description:

Designed for transfusion of blood or blood component
Moulded latex flush ball injection membrane for intermittent medication
Approximately 20 drops / ml
Clear, transparent and soft tubing manufactured from medical grade PVC
With 200 micron filter and 18 G hypodermic needle

Packing:PE package or paper-plastic package

Options available

- With or without Airvent
- With or without Needle
- With or without "Y" site (Latex or Latex - free)
- Luer lock or luer slip

Size:18G/19G/20G/21G/22G/23G/24G/25G/26G/27G

Package:400pcs/ctn, 58*40*38cm



- 1ml (luer)
- 1ml (luer lock)
- 1ml (standard with needle)
- 1ml (long luer lock)
- 1ml (long with needle)

PREFILLABLE SYRINGE

Description:

Quality Standards of glass barrels

- (1) Glass barrel made of neutral-glass clear, type I
- (2) compliance to USP<661> containers-chemical properties- glass containers
- (3) compliance to EP< 3.2.1> glass containers for pharmaceutical use
- (4) compliance to the national glass barrels for prefilled syringes standard.

glass barrels specifications:

1ml long
outside diameter: 8.15 ± 0.1 mm
inside diameter: 6.35 ± 0.1 mm
1ml, 2.25ml, 3ml
outside diameter: 10.85 ± 0.1 mm
inside diameter: 8.65 ± 0.09 mm
Ceramic



- 2.25ml (staked-in needle)
- 2.25ml (cone)
- 2.25ml(luer lock)
- 3ml(luer cone)



PLUNGER ROD

Description:

Prefillable syringe's plunger rod adopt materials with polypropylene, polyethylene, polystyrene

We can produce different colours according with clients requirement



NEEDLE SHIELD

Description:

Quality Standards of Needle Shield And Tip Cap

- (1) Polyisoprene rubber and butyl rubber
- (2) compliance to the national polyisoprene needle shield for prefillable syringes standard.(YBB00102004)
- (3) compliance to USP<381> rubber closures for injection
- (4) compliance to EP< 3.2.9> rubber closures for containers for aqueous parenteral preparations, for powders and for freeze-dried powders

Product Property:

- (1) Well watertightness
- (2) Able to work as sterile protective screening.
- (3) Able to sterilize by ethylene oxide (ethylene oxide gas able to sterilize needlepoint penetrating needle shield)
- (4) Polyisoprene rubber for rigid needle shield and butyl rubber for tip cap



STOPPERS

Description:

Injection need bigger thrust for greater viscous liquid medicine. It's not convenient to hold flange of glass barrels especially to thinner barrel. Backstop developed by Weigao is set on the barrel flange increasing area of the flange, easy to hold syringe, avoid mis-operation and stopper moving before using.



PACKAGING OF PREFILLABLE SYRINGES (1)

Description:

Sterile and ready-to-fill prefillable syringes are supplied to the pharmaceutical industry in special packages. Inside the tubs the syringe barrels are inserted in so-called nests for further processing. This is an ideal system for pharmaceutical production as it is particularly safe and handling is easy.



PACKAGING OF PREFILLABLE SYRINGES (2)

Description:

Sterile and ready-to-fill prefillable syringes are supplied to the pharmaceutical industry in special packages. Inside the tubs the syringe barrels are inserted in so-called nests for further processing. This is an ideal system for pharmaceutical production as it is particularly safe and handling is easy.



STOPPERS PACKAGE (1)


Description:

Quality Standards of rubber plunger

- (1) artificial synthetic butyl rubber
- (2) compliance to the national butyl chloride & butyl bromide rubber plunger for prefillable syringes standard.
- (3) compliance to USP<381> rubber closures for injection
- (4) compliance to EP< 3.2.9> rubber closures for containers for aqueous parenteral preparations, for powders and for freeze-dried powders
- (5) drug industries are responsible for compatibility lab between filling drug and butyl rubber plunger of prefillable syringes.

Product Property:

- (1) compatibility with drug
- (2) easy to push drug (slide well)
- (3) tightness and water tightness are well
- (4) Able to work as sterile protective screening
- (5) Low permeability of gas
- (6) compatible with silicone oil
- (7) easy to insert plunger stopper for filling and packaging machine
- (8) rubber plunger use two kinds of butyl rubber (chlorobutyl , bromobutyl)



STOPPERS PACKAGE (2)

Description:

Quality Standards of rubber plunger

- (1) artificial synthetic butyl rubber
- (2) compliance to the national butyl chloride & butyl bromide rubber plunger for prefillable syringes standard.
- (3) compliance to USP<381> rubber closures for injection
- (4) compliance to EP< 3.2.9> rubber closures for containers for aqueous parenteral preparations, for powders and for freeze-dried powders
- (5) drug industries are responsible for compatibility lab between filling drug and butyl rubber plunger of prefillable syringes.

Product Property:

- (1) compatibility with drug
- (2) easy to push drug (slide well)
- (3) tightness and water tightness are well
- (4) Able to work as sterile protective screening
- (5) Low permeability of gas
- (6) compatible with silicone oil
- (7) easy to insert plunger stopper for filling and packaging machine
- (8) rubber plunger use two kinds of butyl rubber (chlorobutyl , bromobutyl)



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH SMALL WINGS

Description:

F.E.P material
Thin-wall construction permits high flow rates
Radiopaque, x-ray detectable stripes
Clear and transparent
Smooth surfaces
Back-cut special designed needle, siliconized
Rotated wing for fixation and as grip plate
Hydrophobic blood stopper
Color coded
Lock fitting
Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G
Package:30pcs/box,600pcs/ctn,77*53*29cm



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA Y TYPE

Description:

Thin-wall construction permits high flow rates
Radiopaque, x-ray detectable stripes
Clear and transparent
Smooth surfaces
Back-cut special designed needle
Hydrophobic blood stopper
Color coded
Lock fitting
Plastic packing and blister packing are provided
Reduce insertion attempts
Increase protection against blood exposure and needle stick injuries
Simplify IV therapy process
Enhance patient safety and satisfaction

Size: 14G/16G/17G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn,50.5*42.5*28cm



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH INJECTION PORT

Description:

F.E.P material
Thin-wall construction permits high flow rates
Radiopaque, x-ray detectable stripes
Clear and transparent
Smooth surfaces
Back-cut special designed needle, siliconized
Flexible fixation wings
Injection valve with cap
Hydrophobic blood stopper
Color coded
Lock fitting
Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn,48*42*32cm



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA PEN LIKE

Description:

F.E.P material
Thin-wall construction permits high flow rates
Radiopaque, x-ray detectable stripes
Clear and transparent
Smooth surfaces
Back-cut special designed needle, siliconized
Hydrophobic blood stopper
Color coded
Lock fitting
Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn,50*34*33cm

I.V. CANNULA SAFETY TYPE



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

Resoted in the body as combo with infusion sets for transfusion.
This product can be kept in the vein of human for a long time.
It brings less vascular injury. Consequently it largely reduces suffering of patients and reduces intensity of health care work.
The needle retract into the tube after single use, and prevent reusing and protect medical care personnel

Size: 16G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn

I.V. CANNULA SAFETY CLOSED LIGHT PROOF



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

Resoted in the body as combo with infusion sets for transfusion.
This product can be kept in the vein of human for a long time.
It brings less vascular injury. Consequently it largely reduces suffering of patients and reduces intensity of health care work.
The needle retract into the tube after single use, and prevent reusing and protect medical care personnel
Closed light proof to infuse.

I.V. CANNULA PEN-LIKE SAFETY MODEL



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

Thin-wall construction permits high flow rates
Radiopaque,3x-ray detectable stripes
Clear and transparent
Smooth Surfaces
Back-cut special designed needle
Color coded
Lock fitting
Plastic packing and Blister packing are provided
Locking function avoid secondary use,protect the harm of

Size: 16G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn, 48.5*30.5*37.5cm

I.V. CANNULA WITH WINGES SAFETY TYPE



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

Thin-wall construction permits high flow rates
Radiopaque,3x-ray detectable stripes
Clear and transparent
Smooth Surfaces
Back-cut special designed needle
Rotated wing for fixation and as grip plater
Color coded
Lock fitting
Plastic packing and Blister packing are provided
Locking function avoid secondary use,protect the harm of

Size: 16G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn, 60*28.5*52cm

I.V. CANNULA WITH INJECTION VALVE SAFETY TYPE



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

Thin-wall construction permits high flow rates
Radiopaque,3x-ray detectable stripes
Clear and transparent
Smooth Surfaces
Back-cut special designed needle
Flexible fixation wings
Injection valve with cap
Color coded
Lock fitting
Plastic packing and Blister packing are provided
Locking function avoid secondary use,protect the harm of

Size: 16G/18G/20G/22G/24G
Package:50pcs/box,1000pcs/ctn, 41*32*47cm



HEPARIN CAP

Description:

1. Used together with the arterial and venous cannula.
2. Made from fine material (medical grade PVC) and easy to puncture.
3. Latex or latex free
4. It was a tight-fitting adapter.
5. The infusion of Heparin-Sodium can prevent the reflux of blood coagulation
6. Soft colour and very smooth, without any edges and corners.
7. International luer connector, excellent On bio-compatibility.
8. ISO and CE certificate
9. Sterile blister packaging
10. Package: 4000pcs/ctn, 42*26*36cm

HEPARIN CAP WITH CAP

Description:

1. Used together with the arterial and venous cannula.
2. Made from fine material (medical grade PVC) and easy to puncture.
3. Latex or latex free
4. It was a tight-fitting adapter.
5. The infusion of Heparin-Sodium can prevent the reflux of blood coagulation
6. Soft colour and very smooth, without any edges and corners.
7. International luer connector, excellent On bio-compatibility.
8. ISO and CE certificate
9. Sterile blister packaging
10. With cap
11. Package: 4000pcs/ctn, 42*26*36cm





THREE-WAY STOPCOCK

Description:

Rotating two female luer connectors and one male luer.
With or without rotating luer lock.
360°rotation.

Package:1000pcs/ctn, 58*38*38cm



THREE-WAY STOPCOCK WITH EXTENSION TUBE

Description:

Rotating two female luer connectors and one male luer.
With or without rotating luer lock.
360°rotation.
The length of extension tube could be at customer's request.

Three-way stopcock with 10cm extension tube
Three-way stopcock with 15cm extension tube
Three-way stopcock with 20cm extension tube
Three-way stopcock with 25cm extension tube

Package:1000pcs/ctn, 72*26*32cm



TWO-WAY STOPCOCK

Description:

Rotating two female luer connectors and one male luer.
With or without rotating luer lock.
180°rotation.

Package:1000pcs/ctn, 58*38*38cm



TWO-WAY STOPCOCK WITH EXTENSION TUBE

Description:

Rotating two female luer connectors and one male luer.
With or without rotating luer lock.
180°rotation.
The length of extension tube could be at customer's request.

Two-way stopcock with 10cm extension tube
Two-way stopcock with 15cm extension tube
Two-way stopcock with 20cm extension tube
Two-way stopcock with 25cm extension tube

Package:1000pcs/ctn, 72*26*32cm



LARYNGEAL MASK Silicone Disposable

Description:

- Disposable Laryngeal Mask Airway.
- 100% medical-grade Dow Corning silicone.
- Three colors for the cuff: light blue, light yellow, clear.
- Smooth connecting, very high quality.
- Package: 2pcs/box, 50pcs/carton, 77*40*27 cm

Disposable yellow cuff with aperture bars	Disposable blue cuff with aperture bars	Disposable clear cuff with aperture bars	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
AP10-Y1	AP10-B1	AP10-C	1	4	0-5
AP15-Y1	AP15-B1	AP15-C	1.5	7	5-10
AP20-Y1	AP20-B1	AP20-C	2	10	10-20
AP25-Y1	AP25-B1	AP25-C	2.5	14	20-30
AP30-Y1	AP30-B1	AP30-C	3	20	30-50
AP40-Y1	AP40-B1	AP40-C	4	30	50-70
AP50-Y1	AP50-B1	AP50-C	5	40	70-100



LARYNGEAL MASK Silicone Reusable

Description:

- Reusable Laryngeal Mask Airway.
- 100% medical-grade Dow Corning silicone.
- Both with and without aperture bars are available.
- Can be autoclaved repeatedly at 134°C (deflate the cuff before sterilization).
- Reusable up to 40 times.
- Package: 2pcs/box, 50pcs/carton, 77*40*27 cm

Reusable yellow cuff with aperture bars	Reusable blue cuff with aperture bars	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
AP10-Y2	AP10-B2	1	4	0-5
AP15-Y2	AP15-B2	1.5	7	5-10
AP20-Y2	AP20-B2	2	10	10-20
AP25-Y2	AP25-B2	2.5	14	20-30
AP30-Y2	AP30-B2	3	20	30-50
AP40-Y2	AP40-B2	4	30	50-70
AP50-Y2	AP50-B2	5	40	70-100



LARYNGEAL MASK Reinforced Silicone Disposable

Description:

- Reinforced Laryngeal Mask Airway, disposable.
- Made of 100% medical-grade Dow Corning silicone.
- A stainless steel spring is inserted in the tube to reduce the risk of kinking.
- Three colors for the cuff: light blue, light yellow, clear.
- Package: 2pcs/box, 50pcs/carton, 77*40*27 cm

Silicone disposable reinforced clear cuff with aperture bars	Silicone disposable reinforced clear cuff without aperture bars	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
AP10-C2	AP10-C3	1	4	0-5
AP15-C2	AP15-C3	1.5	7	5-10
AP20-C2	AP20-C3	2	10	10-20
AP25-C2	AP25-C3	2.5	14	20-30
AP30-C2	AP30-C3	3	20	30-50
AP40-C2	AP40-C3	4	30	50-70
AP50-C2	AP50-C3	5	40	70-100



LARYNGEAL MASK PVC Disposable

Description:

- PVC disposable Laryngeal Mask Airway.
- Made of medical-grade PVC.
- Color - coded, easy to identify sizes.

Package:

- **Paper-poly Pouch:**
2pcs/box, 50pcs/carton, 77*40*27cm
- **Blister Pack:**
10pcs/box, 50pcs/carton
1#~2.5#: 47.5*33*28cm
3.0#, 4.0#: 51.5*41*33cm
5#: 55*43.5*37cm

Type	Product No	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
Standard	AP10-1	1	4	0-5
	AP15-1	1.5	7	5-10
	AP20-1	2	10	10-20
	AP25-1	2.5	14	20-30
	AP30-1	3	20	30-50
	AP40-1	4	30	50-70
	AP50-1	5	40	70-100



LARYNGEAL MASK Reinforced PVC Disposable

Description:

- Reinforced PVC Laryngeal Mask Airway.
- Made of medical-grade PVC.
- A stainless steel spring is inserted in the tube to reduce the risk of kinking.
- Color - coded, easy to identify sizes.

Package:


- **Paper-poly Pouch:**
2pcs/box, 50pcs/carton, 77*40*27cm
- **Blister Pack:**
10pcs/box, 50pcs/carton
1#~2.5#: 47.5*33*28cm
3.0#, 4.0#: 51.5*41*33cm
5#: 55*43.5*37cm

Type	Product No	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
Reinforced	AP10-1	1	4	0-5
	AP15-1	1.5	7	5-10
	AP20-1	2	10	10-20
	AP25-1	2.5	14	20-30
	AP30-1	3	20	30-50
	AP40-1	4	30	50-70
	AP50-1	5	40	70-100

ENDOTRACHEAL TUBE Oral/Nasal

Description:

- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- Murphy eye to avoid complete respiratory obstruction.
- Radiopaque line throughout the tube for X-ray visualization.
- Package:20pcs/box, 200pcs/carton, 60*37*34cm



Type	Product No	Size(#)
Oral/Nasal Without Cuff	AP20C	2
	AP25C	2.5
	AP30C	3
	AP35C	3.5
	AP40C	4
	AP45C	4.5
	AP50C	5
	AP55C	5.5
	AP60C	6
	AP65C	6.5
	AP70C	7
	AP75C	7.5
	AP80C	8
	AP85C	8.5
AP90C	9	
AP95C	9.5	
AP10C	10	

Type	Product No	Size(#)
Oral/Nasal With Cuff	AP30C	3
	AP35C	3.5
	AP40C	4
	AP45C	4.5
	AP50C	5
	AP55C	5.5
	AP60C	6
	AP65C	6.5
	AP70C	7
	AP75C	7.5
	AP80C	8
	AP85C	8.5
	AP90C	9
	AP95C	9.5
AP10C	10	

REINFORCED ENDOTRACHEAL TUBE Oral/Nasal

Description:

- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- Murphy eye to avoid complete respiratory obstruction.
- Radiopaque line throughout the tube for X-ray visualization.
- A stainless spring is inserted in the tube to minimize the risk of kinking or crushing.
- Straight reinforced endotracheal tube with preloaded stylet is very convenient for use.
- Package:20pcs/box, 200pcs/carton, 60*37*34cm



Type	Reinforced with cuff-curved	Reinforced without cuff-curved	Reinforced with cuff-straight with preloaded stylet	Reinforced without cuff-straight with preloaded stylet	Size(#)
Both curved and straight reinforced endotracheal tubes are available. Preloaded stylet is optional.	AP30-1	AP30-2	AP30-3	AP30-4	3
	AP35-1	AP35-2	AP35-3	AP35-4	3.5
	AP40-1	AP40-2	AP40-3	AP40-4	4
	AP45-1	AP45-2	AP45-3	AP45-4	4.5
	AP50-1	AP50-2	AP50-3	AP50-4	5
	AP55-1	AP55-2	AP55-3	AP55-4	5.5
	AP60-1	AP60-2	AP60-3	AP60-4	6
	AP65-1	AP65-2	AP65-3	AP65-4	6.5
	AP70-1	AP70-2	AP70-3	AP70-4	7
	AP75-1	AP75-2	AP75-3	AP75-4	7.5
	AP80-1	AP80-2	AP80-3	AP80-4	8
	AP85-1	AP85-2	AP85-3	AP85-4	8.5
	AP90-1	AP90-2	AP90-3	AP90-4	9
	AP95-1	AP95-2	AP95-3	AP95-4	9.5
	AP10-1	AP10-2	AP10-3	AP10-4	10



NASAL PREFORMED TUBE

Description:

Available with cuff and without cuff, and reinforced tube available

Made of medical grade transparent thermosensitive PVC.
Cuffed tube with high volume low pressure cuff.
Special formed for nasal use.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package:20pcs/ctn, 200pcs/ctn, 60*37*34cm



ORAL PREFORMED TUBE

Description:

Available with cuff and without cuff, and reinforced tube available

Made of medical grade transparent thermosensitive PVC.
Cuffed tube with high volume low pressure cuff.
Special formed for oral use.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package:20pcs/box, 200pcs/ctn, 60*37*34cm



INJECTABLE ENDOTRACHEAL TUBE

Description:

Made of medical grade PVC and latex free.

With injection port

Size: 4.5/5.0/5.5/6.0/6.5/7.0/7.5/8.0/8.5/9.0

Package:10pcs/box,100pcs/ctn,45*39.5*32cm



ENDOBONCHIAL TUBE

Description:

Three styles of bronchial cuffs are available.


Two styles of connectors, fixed and non-fixed.

Low-pressure cuffs can help to reduce the lesion of mucosa.

Available also in a set together with connector and three pieces of suction catheters.

Size: Right 26Fr/28Fr/32Fr/35Fr/37Fr/39Fr/41Fr
Left 26Fr/28Fr/32Fr/35Fr/37Fr/39Fr/41Fr

Package:10pcs/ctn,53*33*30cm



TRACHEOTOMY TUBE

Description:

Formed in 100±10° C.

With two pieces of 65cm ribbons.

Radiopaque line throughout the tube.

Made of transparent, nontoxic PVC and latex free.

Available with cuff(high volume low pressure)and without cuff.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package:10pcs/box,100pcs/ctn,45*39.5*32cm



REINFORCED TRACHEOTOMY TUBE

Description:

Made of medical grade PVC and latex free.

Spiral reinforcement minimizes crushing or kinking.

Available with cuff(high volume low pressure)and without cuff.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package:10pcs/box,100pcs/ctn,45*39.5*32cm



GUEDEL AIRWAY - with colored mark

Description:

- Color - coded Guedel Airway, for easy identification of sizes.
- Made of polypropylene.
- Package:500pcs/carton, 62*47*32cm

Product No	Size(#)
AP40	40MM
AP50	50MM
AP60	60MM
AP70	70MM
AP80	80MM
AP90	90MM
AP100	100MM
AP110	110MM
AP120	120MM



ENDOTRACHEAL TUBE HOLDER

Description:

- This Endotracheal Tube Holder fits any size of endotracheal tubes.
- Patent product.
- Package:100pcs/carton, 48.5*33*41cm

Product No	Size(#)
AP40	Adult
AP50	Pediatric



BREATHING CIRCUIT - CORRUGATED

Description:

- Breathing Circuit, corrugated tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
- Size(#): Adult/Pediatric
- Remark:Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
- Package: 500pcs/carton, 65*46*35cm



BREATHING CIRCUIT - COAXIAL

Description:

- Coaxial breathing circuit, standard connectors (15mm, 22mm).
- Inside gas sampling line. Very high quality.
- The tube color can be customized.
- Size(#): Adult/Pediatric
- Remark:Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
- Package:500pcs/carton, 65*46*35cm



BREATHING CIRCUIT - SMOOTHBORE

Description:

- Breathing Circuit, smoothbore tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
- Size(#): Adult/Pediatric
- Remark:Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
- Package:500pcs/carton, 65*46*35cm



BREATHING CIRCUIT - EXPANDABLE

Description:

- Breathing Circuit, expandable tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
- Size(#): Adult/Pediatric
- Remark:Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
- Package:500pcs/carton, 65*46*35cm



ANESTHESIA MASK

Description:

- Anesthesia Mask, made of clear, soft, medical-grade PVC.
- Seal tight against mouth.
- Color coded to easily identify sizes.
- Size(#): 1, 2, 3, 4, 5, 6
- Package: 100pcs/carton, 55*37*34cm



BREATHING CIRCUIT - CUSTOMIZED

Description:

- Our breathing circuit (ventilator circuit, anesthesia circuit) can be customized at many lengths, and be equipped with Water Trap, Breathing Bag (latex or latex free), or Filter, HMEF, Catheter Mount, etc. If you have any specific requirements, please just contact us.
- We have Corrugated, Expandable, Smoothbore, Coaxial and other tubes, various types.



BREATHING BAG

Description:

- Size: 0.5L, 1L, 2L, 3L
- Equipped with 22F standard connector.
- Both Latex and Latex-free are available.
- Color: grey and blue are latex, and green is latex-free.



NASAL OXYGEN CANNULA

Description:

- Nasal Cannula, made of clear, medical-grade PVC.
- Diameter 6mm, kinking resistant, high quality.
- Size(#):

Adult straight	Adult soft curved tip
Child straight	Child soft curved tip
Infant straight	Infant soft curved tip
Neonate straight	Neonate soft curved tip
- Package: 100pcs/carton, 35*37*34cm



OXYGEN MASK


Description:

Made of medical PVC.
Adjustable nose clip.
Medium concentration.
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

C

Package:100pcs/ctn, 50*37*34cm



NEBULIZER MASK

Description:

Made of medical PVC.
Without nose clip.
Medium concentration.
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

Package:100pcs/ctn, 50*37*34cm



OXYGEN MASK WITH BAG

Description:

Made of medical PVC.
Adjustable nose clip.
High concentration.
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

Package:100pcs/ctn, 50*37*34cm



VENTURI MASK

Description:

Made of medical PVC.
Adjustable nose clip.
Oxygen concentration is 24%-50%.
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

Package:100pcs/ctn, 50*37*34cm




ADJUSTABLE VENTURI MASK

Description:

Made of medical PVC.
Adjustable nose clip.
Venturi valve:blue 24%,yellow 28%,white 31%,green 35%,pink 40%,orange 50%.
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

Package:100pcs/ctn, 50*37*34cm




SUCTION CATHETER

Description:

- Made from clear soft PVC.
- Four kind of connector available with T-type connector, Y-type connector, Cap-cone connector, Y-type transparent connector.
- Colour coded for size identification.
- Can With Scale line.
- Length: 500mm.

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F221 05000	5	Light Gray
APF-2	F221 06000	6	Light Green
APF-3	F250 08000	8	Blue
APF-4	F250 10000	10	Black
APF-5	F250 12000	12	White
APF-6	F250 14000	14	Green
APF-7	F250 16000	16	Orange
APF-8	F250 18000	18	Red
APF-9	F250 20000	20	Yellow
APF-10	F250 22000	22	Blue



RECTAL TUBE

Description:

- Soft, frosted and kink resistant PVC tubing.
- For introduction of enema solution into rectum to release
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage
- Colour coded plain connector for easy identification of size
- Length:400mm

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F260 2400	24	Violet
APF-2	F260 2600	26	Red
APF-3	F260 2800	28	Yellow
APF-4	F260 3000	30	Black
APF-5	F260 3200	32	Blue
APF-6	F260 3400	34	Green
APF-7	F260 3600	36	Light Blue



Female Type-200mm



Standard-400mm

NELATON CATHETER -FEMALE TYPE/STANDARD

Description:

- Designed for short term and temporary drainage of urine.
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage
- Legth:400mm/200mm
- Package: 55*27*37cm
- Female Type Package : 1000pcs, 55*29*37cm

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code	Pc/Ctn
APF-1	F250 0600	6	Light Green	800pcs
APF-2	F250 0800	8	Blue	800pcs
APF-3	F250 1000	10	Black	600pcs
APF-4	F250 1200	12	White	600pcs
APF-5	F250 1400	14	Green	600pcs
APF-6	F250 1600	16	Orange	500pcs
APF-7	F250 1800	18	Red	500pcs
APF-8	F250 2000	20	Yellow	500pcs
APF-9	F250 2200	22	Blue	500pcs
APF-10	F250 2400	24	Violet	500pcs



YANKAUER SUCTION HANDLE

Description:

- Made of K-resin or PMMA
- Crown tip & standard tip
- With or without vent
- Available with suction connection tube



SUCTION CONNECTING TUBE

Description:

- Made of non-toxic PVC, medical grade
- Hex-arris designed to avoid blocking due to high pressure
- Suitable for pre-operative removal of secretion and body fluids
- Standard connector at distal end of tubing to connect Yankaur Handle
- Tubing measures 1/4 or 3/16 in diameter, different tube length available



THORACIC DRAINAGE TUBE

Description:

- Non-toxic, medical grade PVC material
- X-ray opaque line.
- Type: round or fish-mouth tip
- Length: 500mm
- Size: 20/22/24/26/28/30/32/34/36



ABDOMINAL DRAINAGE TUBE

Description:

- Non-toxic, medical grade PVC material
- X-ray opaque line.
- Length: 500mm
- Size: 20/22/24/26/28/30/32/34/36



FEEDING TUBE MULTIFUNCTION CONNECTOR/LUER FITTING

Description:

- Material: Made form Medical grade PVC
- Two type: standard connector and Multifunction connector
- Atraumatic, soft, rounded, opened tip with two lateral eyes .
- With X-ray line available
- F4~F8 Length:40cm; F10~F22 Length:120cm

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F221 04000	4	Red
APF-2	F221 05000	5	Light Gray
APF-3	F221 06000	6	Light Green
APF-4	F250 08000	8	Blue
APF-5	F250 10000	10	Black
APF-6	F250 12000	12	White
APF-7	F250 14000	14	Green
APF-8	F250 16000	16	Orange
APF-9	F250 18000	18	Red
APF-10	F250 20000	20	Yellow
APF-11	F250 22000	22	Blue
APF-12	F250 24000	24	Violet



Ryle's Stomach Tube



Standard Type

STOMACHE TUBE RYLE'S STOMACH TUBE /STANDARD TYPE


Description:

- X- ray line/ radio-opaque line is available for accurate placement
- Tube end is sealed with stainless steel balls, to assist the catheter indwelling in a certain place

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F250 1200	12	White
APF-2	F250 1400	14	Green
APF-3	F250 1600	16	Orange
APF-4	F250 1800	18	Red
APF-5	F250 2000	20	Yellow
APF-6	F250 2200	22	Blue
APF-7	F250 2400	24	Violet



2-way



3-way

SILICON COATED LATEX FOLEY CATHETERS

Description:

Made of from natural latex, silicone coated
Rubber valve and plastic valve for different customer needs
Length 400mm
Package: 500pcs/ctn, 52.5x41x43cm

2-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP2FC-1	12	White	5-15	Without stylet
AP2FC-2	14	Green		
AP2FC-3	16	Orange		
AP2FC-4	18	Red		
AP2FC-5	20	Yellow		
AP2FC-6	22	Purple		
AP2FC-7	24	Blue		
AP2FC-8	26	Pink		
AP2FC-9	16	Orange	30	
AP2FC-10	18	Red		
AP2FC-11	20	Yellow		
AP2FC-12	22	Purple		
AP2FC-13	24	Blue		
AP2FC-14	26	Pink		

3-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP3FC-1	16	Orange	30	3-way
AP3FC-2	18	Red		
AP3FC-3	20	Yellow		
AP3FC-4	22	Purple		
AP3FC-5	24	Blue		
AP3FC-6	26	Pink		

ALL SILICON FOLEY CATHETERS

Description:

Made of medical-grade silicone
Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently
Length 400mm
Package: 500pcs/ctn, 52x43.5x39.5cm

2-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP2FC-1	12	White	5-15	Without stylet
AP2FC-2	14	Green		
AP2FC-3	16	Orange		
AP2FC-4	18	Red		
AP2FC-5	20	Yellow		
AP2FC-6	22	Purple		
AP2FC-7	24	Blue		
AP2FC-8	26	Pink		
AP2FC-9	16	Orange	30	
AP2FC-10	18	Red		
AP2FC-11	20	Yellow		
AP2FC-12	22	Purple		
AP2FC-13	24	Blue		
AP2FC-14	26	Pink		

3-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP3FC-1	16	Orange	30	3-way
AP3FC-2	18	Red		
AP3FC-3	20	Yellow		
AP3FC-4	22	Purple		
AP3FC-5	24	Blue		






PEDIATRIC URINE BAG

Description:


Size: 100ml/150ml/200ml
Safety Adhesive, Sterile, Latex Free



URINE BAG WITHOUT OUTLET

Description:

Size: 2000ml
Packing:
1pcs/polybag
10pcs/big PE bag
250pcs/ctn
54x40x27cm



URINE BAG WITH PUSHING VALVE

Description:

Size: 2000ml
Packing :
1pcs/polybag
10pcs/big PE bag
250pcs/ctn
54x40x27cm



URINE BAG WITH SCREW VALVE

Description:

Size: 2000ml
Packing:
1pcs/polybag
10pcs/big PE bag
250pcs/ctn
54x40x27cm



URINE BAG WITH T VALVE

Description:

Size: 2000ml
Packing :
1pcs/polybag
10pcs/big PE bag
250pcs/ctn
54x40x27c



(1)



(2)



(3)



(4)

TWIST TOP BLOOD LANCET

Description:

- Be sterilized by GAMMA
- Blood lancet is a small medical implement used for capillary blood sampling.
- Blood Lancets are used to make punctures, such as a finger stick, to obtain small blood specimens. Blood lancets are generally disposable.
- Blood Lancets are also used to prick the skin in skin testing for allergies.
- Blood lancet is used with Lancing Device.
- We can process the products according to the samples provided by our customers.

(1) Size:21G,23G,26G,28G,30G

Package: 100pcs/box
200box/ctn
49*35*35cm

(2) Size:28G,30G

Package: 200pcs/box
100box/ctn
40*27*37cm

(3) Size:28G,30G

Package: 200pcs/box
100box/ctn
46*42.5*29cm

(4) Size:28G,30G

Package: 200pcs/box
100box/ctn
46*42.5*29 cm



(1)



(2)

LANCING DEVICES

Description:

A blood-sampling device, also known as a lancing device, is a reusable instrument equipped with a lancet. It is also most commonly used by diabetics during blood glucose monitoring. The depth of skin penetration can be adjusted for various skin thicknesses. Long lancing devices are used for fetal scalp blood testing to get a measure of the acid base status of the fetus.

(1) Package: 1pcs/box, 100box/ctn , 27.5*22*22cm

(2) Package: 1pcs/box, 100box/ctn , 27*26.5*14.5cm



1 Quick, easy lancet insertion and disposal.



2 It is simple to remove and twist off lancet protective cap.



3 Five depth settings for optimum skin penetration.



4 Sliding barrel signals readiness with a soft click.



5 One-handed lancing with comfortable trigger button.

STAINLESS STEEL BLOOD LANCET

Description:

- Be sterilized by GAMMA
- Blood lancet is a small medical implement used for capillary blood sampling.
- Blood Lancets are used to make punctures, such as a finger stick, to obtain small blood specimens. Blood lancets are generally disposable.
- Blood Lancets are also used to prick the skin in skin testing for allergies.
- We can process the products according to the samples provided by our customers.

Size: L, S

Package: 200pcs/box, 100box/ctn , 33.5*29*29.5cm



SAFE BLOOD LANCET

Description:

Safe Lancets can maximize avoid the possibility of health care workers worldwide operation in blood samples infected, such as HIV, hepatitis and other diseases. After the launch lancets, needle core will be locked, so that the lancet can only be used once, you can ensure the safety of users; High-quality, ultra-sharp edged needle head design that can quickly pierce the skin, reducing the patient's pain; A variety of models and specifications of needles piercing depth, suitable for the vast majority of blood collection absolutely need; Suitable for diabetic blood sugar and other tests and diagnostics

(1) Size:21G,23G,26G(2.4MM)


Package: 100pcs/box, 20box/ctn , 39.5*22*26cm



(1)



1 Push the lancet cap into the body of lancet until it clicks.



2 Carefully twist off the protective cap until it is separated from the device.



3 Push the lancet firmly onto the chosen site, then press trigger.




4 Dispose of the lancet in a suitable sharps container.




(2)

(2) Size:21G(2.4MM) 23G(2.4MM) 26G(1.8MM,2.4MM)


Package: 100pcs/box, 20box/ctn , 52.5*32.5*17cm



1 Carefully pull off the protective cap.



2 Push the lancet firmly onto the chosen site to fire.



3 Dispose of lancet in a suitable sharps container.



DISPOSABLE VACUUM BLOOD COLLECTION TUBES

Description:

Material PET/PP/GLASS
Cap color Blue
Additive 3.2% Sodium Citrate / 3.8% Sodium Citrate(1:9)
Specimen Whole Blood or plasma
Specs 13*75mm
Invert the tube for 6-8 times gently.
Size:1.8ml/2.7ml/3.6ml/4.5ml
Package:100pcs*12racks/carton, 51*36*21cm



DISPOSABLE VACUUM BLOOD COLLECTION TUBES

Description:

Material PET/PP/GLASS
Cap color Purple
Additive EDTA K2 / EDTA K3 / EDTA 2Na
Specimen Whole Blood
Specs 13*75mm
Invert the tube for 6-8 times immediately after drawing blood. When the outside temperature is higher than 35°C, the inverting times should be increased accordingly.
Size:1.0ml/2.0ml/3.0ml/4.0ml /5.0ml/6.0ml/7.0ml
Package:100pcs*12racks/carton, 51*36*21cm



DISPOSABLE VACUUM BLOOD COLLECTION TUBES

Description:

Material PET/PP/GLASS
Cap color Black
Additive 3.8% Sodium Citrate(1:4)
Specimen Whole Blood
Specs 120*8mm
Invert the tube for 6-8 times gently.
Size:1.28ml
Package:100pcs*12racks/carton, 51*36*21cm



DISPOSABLE VACUUM BLOOD COLLECTION TUBES

Description:

Material PET/PP/GLASS
Cap color Green
Additive Lithium Heparin / Sodium Heparin
Specimen Plasma
Specs 13*75mm
Draw required volume and invert the tube for 6-8 times gently.
Size:1.0ml/2.0ml/3.0ml/4.0ml /5.0ml/6.0ml/7.0ml
Package:100pcs*12racks/carton, 51*36*21cm



DISPOSABLE VACUUM BLOOD COLLECTION TUBES

Description:

Material PET/PP/GLASS
Cap color Yellow
Additive Gel & Clot Activator
Specimen Serum
Specs 13*75mm,13*100mm,16*100mm
Draw volume 1-9ml
The Gel is infused in the middle of the tube, and a hole penetrates through the gel vertically. The new technics will reduce the time of centrifuge and the relative centrifugal force, meanwhile, high quality serum can be separated well.
Size:1.0ml/2.0ml/3.0ml/4.0ml /5.0ml/6.0ml/7.0ml/8.0ml/9.0ml
Package:100pcs*12racks/carton, 51*36*21cm



PLASTIC BLOOD BAG FOR SINGLE USE

Description:

- 1.Primary Bag capacity available in 150ml to 500ml
- 2.Satellite Bag capacity available in 100ml to 450ml
- 3.Package: Each set packed in one PE or PET compounded vacuum pouch, 3 to 5 sets packed in one aluminum foil compounded pouch, and then 30 to 50 sets packed in one carton

OPTIONAL FEATURES

1. Attached Blood Sampling Bag
2. Needle Protector
3. Tube Holder
4. "Y" injection port
5. In-line Leukocyte Reduction Filter
6. Platelet Storage Bag (platelet storage for 3 to 5 days or up to 5 days)
7. Type of outlet Ports (twist off or tear off)
8. Volume of primary bag and satellite bags
9. Rolled Bag or Blow-extruded Bag

Size: Needle 16G/17G

Measure of carton :

- One 250ml Bag with CPDA1 and one 100 to 200ml empty satellite bag
- One 350ml Bag with CPDA1 and one 200 to 300ml empty satellite bag
- One 450ml Bag with CPDA1 and one 300 to 400ml empty satellite bag
- One 500ml Bag with CPDA1 and one 300 to 400ml empty satellite bag



BLOOD COLLECTION NEEDLE

Description:

Used to collect blood for testing, usually used together with vacuum tube
It provide fast and convenient blood multiple sampling by single puncture.
High quality ultra-sharp needle assure atraumatic insertion and painless penetration
Color-coded cap and hub assure accurate and convenient gauge size indentification
Safety sleeve for this model is available to protect sampling person from accidental hurt or cross-infection
Size: 20G,21G,22G,23G,25G



BLOOD COLLECTION NEEDLE

Description:

Used to collect blood for testing, usually used together with vacuum tube
It provide fast and convenient blood multiple sampling by single puncture.
High quality ultra-sharp needle assure atraumatic insertion and painless penetration
Color-coded cap and hub assure accurate and convenient gauge size indentification
Safety sleeve for this model is available to protect sampling person from accidental hurt or cross-infection
Size: 20G,21G,22G,23G,25G