



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



LARYNGEAL MASK Silicone Disposable

Description:

- Disposable Laryngeal Mask Airway.
- 100% medical-grade Dow Corning silicone.
- Three colors for the cuff: light blue, light yellow, clear.
- Smooth connecting, very high quality.
- Package: 2pcs/box, 50pcs/carton, 77*40*27 cm

Disposable yellow cuff with aperture bars	Disposable blue cuff with aperture bars	Disposable clear cuff with aperture bars	Size(#)	Maximum Cuff Inflation Volumes (ml)	Patient Weight (KG)
AP10-Y1	AP10-B1	AP10-C	1	4	0-5
AP15-Y1	AP15-B1	AP15-C	1.5	7	5-10
AP20-Y1	AP20-B1	AP20-C	2	10	10-20
AP25-Y1	AP25-B1	AP25-C	2.5	14	20-30
AP30-Y1	AP30-B1	AP30-C	3	20	30-50
AP40-Y1	AP40-B1	AP40-C	4	30	50-70
AP50-Y1	AP50-B1	AP50-C	5	40	70-100



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



LARYNGEAL MASK Silicone Reusable

Description:

- Reusable Laryngeal Mask Airway.
- 100% medical-grade Dow Corning silicone.
- Both with and without aperture bars are available.
- Can be autoclaved repeatedly at 134°C (deflate the cuff before sterilization).
- Reusable up to 40 times.
- Package: 2pcs/box, 50pcs/carton, 77*40*27cm

Reusable yellow cuff with aperture bars	Reusable blue cuff with aperture bars	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
AP10-Y2	AP10-B2	1	4	0-5
AP15-Y2	AP15-B2	1.5	7	5-10
AP20-Y2	AP20-B2	2	10	10-20
AP25-Y2	AP25-B2	2.5	14	20-30
AP30-Y2	AP30-B2	3	20	30-50
AP40-Y2	AP40-B2	4	30	50-70
AP50-Y2	AP50-B2	5	40	70-100



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



LARYNGEAL MASK Reinforced Silicone Disposable

Description:

- Reinforced Laryngeal Mask Airway, disposable.
- Made of 100% medical-grade Dow Corning silicone.
- A stainless steel spring is inserted in the tube to reduce the risk of kinking.
- Three colors for the cuff: light blue, light yellow, clear.
- Package: 2pcs/box, 50pcs/carton, 77*40*27cm

Silicone disposable reinforced clear cuff with aperture bars	Silicone disposable reinforced clear cuff without aperture bars	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
AP10-C2	AP10-C3	1	4	0-5
AP15-C2	AP15-C3	1.5	7	5-10
AP20-C2	AP20-C3	2	10	10-20
AP25-C2	AP25-C3	2.5	14	20-30
AP30-C2	AP30-C3	3	20	30-50
AP40-C2	AP40-C3	4	30	50-70
AP50-C2	AP50-C3	5	40	70-100



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



LARYNGEAL MASK PVC Disposable

Description:

- PVC disposable Laryngeal Mask Airway.
- Made of medical-grade PVC.
- Color - coded, easy to identify sizes.

Package:

- **Paper-poly Pouch:**

2pcs/box, 50pcs/carton, 77*40*27cm

- **Blister Pack:**

10pcs/box, 50pcs/carton

1#~2.5#: 47.5*33*28cm

3.0#, 4.0#: 51.5*41*33cm

5#: 55*43.5*37cm

Type	Product No	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
Standard	AP10-1	1	4	0-5
	AP15-1	1.5	7	5-10
	AP20-1	2	10	10-20
	AP25-1	2.5	14	20-30
	AP30-1	3	20	30-50
	AP40-1	4	30	50-70
	AP50-1	5	40	70-100



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



LARYNGEAL MASK Reinforced PVC Disposable

Description:

- Reinforced PVC Laryngeal Mask Airway.
- Made of medical-grade PVC.
- A stainless steel spring is inserted in the tube to reduce the risk of kinking.
- Color - coded, easy to identify sizes.

Package:

- **Paper-poly Pouch:**
2pcs/box, 50pcs/carton, 77*40*27cm
- **Blister Pack:**
10pcs/box, 50pcs/carton
1#~2.5#: 47.5*33*28cm
3.0#, 4.0#: 51.5*41*33cm
5#: 55*43.5*37cm

Type	Product No	Size(#)	Maximum Cuff Inflation Volumes(ml)	Patient Weight (KG)
Reinforced	AP10-1	1	4	0-5
	AP15-1	1.5	7	5-10
	AP20-1	2	10	10-20
	AP25-1	2.5	14	20-30
	AP30-1	3	20	30-50
	AP40-1	4	30	50-70
	AP50-1	5	40	70-100