



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



SPINAL NEEDLES COMMON POINT

Description:

Latex Free
Convenient blister packaging
Siliconized for smooth comfortable
Penetration
Transparent plastic hub

Package: 50pcs/box, 20box/carton, 1000pcs/ctn,
51cm*42cm*25.5cm

Description

16G(1.60mm)X3 1/2"
17G(1.40mm)X3 1/2"
18G(1.20mm)X3 1/2"
19G(1.10mm)X3 1/2"
20G(0.90mm)X3 1/2"
21G(0.80mm)X3 1/2"
22G(0.70mm)X3 1/2"
23G(0.60mm)X3 1/2"
24G(0.55mm)X3 1/2"
25G(0.50mm)X3 1/2"
26G(0.45mm)X3 1/2"
27G(0.40mm)X3 1/2"

Hub Color

□ White
■ Purple
■ Pink
■ Cream
■ Yellow
■ Deep Green
■ Black
■ Deep Blue
■ Violet
■ Orange
■ Brown
■ Grey



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



SPINAL NEEDLES PENCIL POINT

Description:

Latex Free
Convenient blister packaging
Siliconized for smooth comfortable
Penetration
Transparent plastic hub

Package: 50pcs/box, 20box/carton, 1000pcs/ctn,
51cm*42cm*25.5cm

Description

16G(1.60mm)X3 1/2"
17G(1.40mm)X3 1/2"
18G(1.20mm)X3 1/2"
19G(1.10mm)X3 1/2"
20G(0.90mm)X3 1/2"
21G(0.80mm)X3 1/2"
22G(0.70mm)X3 1/2"
23G(0.60mm)X3 1/2"
24G(0.55mm)X3 1/2"
25G(0.50mm)X3 1/2"
26G(0.45mm)X3 1/2"
27G(0.40mm)X3 1/2"

Hub Color

White
 Purple
 Pink
 Cream
 Yellow
 Deep Green
 Black
 Deep Blue
 Violet
 Orange
 Brown
 Grey