



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



NASOPHARYNGEAL AIRWAY

Description:

Slidable rings are design to adjust depth of the tube for more comfort while also prevent tube from slipping the nasal airway beyond the desired depth.

Thickened flange and trumpet flange are designed to prevent tube slipping into the nasal airway.

Cat. No.	Size ID (mm)	Cat. No.	Size ID (mm)
29100025	2.5	29100065	6.5
29100030	3.0	29100070	7.0
29100035	3.5	29100075	7.5
29100040	4.0	29100080	8.0
29100045	4.5	29100085	8.5
29100050	5.0	29100090	9.0
29100055	5.5	29100095	9.5
29100060	6.0	29100100	10.0