




青岛爱普国际贸易有限公司 ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-532-5568 9697 Fax:86-532-5568 9611 E-mail:interventional@allpro.com.cn www.allprocorporation.com



PTCA Balloon Dilatation Catheter CE 0197



- 0.017" small tip profile for excellent crossability.
- Soft tapered tip, provides outstanding ability to encircled the guidewire, minimizing vessel injury.
Stable optimized lubricious hydrophilic coating for smooth delivery.
- Transparent hub is convenient to observe the air bubble, reducing the operation risk.
- The excellent smooth transition between hub and proximal shaft is designed for enhanced pushability and torque transfer.

Technical Parameters

Catheter Design	Rapid Exchange		
Balloon Compliance	Semi-compliant		
Marker	2		
Material	Nylon 12		
Nominal Pressure	1.5~3.0mm 8atm	3.5~4.0mm	6atm
Rated Burst Pressure	1.5~3.5mm 16atm	4.0mm	14atm
Distal Shaft (outer diameter)	1.5~2.0mm 2.6F	2.5~4.0mm	2.8F
Proximal Shaft (outer diameter)	1.5~2.0mm 2.1F	2.5~4.0mm	2.3F
Crossing Profile	0.026~0.030"		
Usable Catheter Length	135cm		
Minimum Guiding Catheter Compatibility	5Fr/0.056"		
Guidewire Compatible	0.014"		

Ordering Information

Length (mm)	Diameter(mm)							
	1.5	1.75	2.0	2.5	2.75	3.0	3.5	4.0
9	PRX15009	PRX17009	PRX20009	PRX25009	PRX27009	PRX30009	PRX35009	PRX40009
11	PRX15011	PRX17011	PRX20011	PRX25011	PRX27011	PRX30011	PRX35011	PRX40011
12	PRX15012	PRX17012	PRX20012	PRX25012	PRX27012	PRX30012	PRX35012	PRX40012
14	PRX15014	PRX17014	PRX20014	PRX25014	PRX27014	PRX30014	PRX35014	PRX40014
16	PRX15016	PRX17016	PRX20016	PRX25016	PRX27016	PRX30016	PRX35016	PRX40016
17	PRX15017	PRX17017	PRX20017	PRX25017	PRX27017	PRX30017	PRX35017	PRX40017
18	PRX15018	PRX17018	PRX20018	PRX25018	PRX27018	PRX30018	PRX35018	PRX40018
20	PRX15020	PRX17020	PRX20020	PRX25020	PRX27020	PRX30020	PRX35020	PRX40020
22	PRX15022	PRX17022	PRX20022	PRX25022	PRX27022	PRX30022	PRX35022	PRX40022
23	PRX15023	PRX17023	PRX20023	PRX25023	PRX27023	PRX30023	PRX35023	PRX40023
24	PRX15024	PRX17024	PRX20024	PRX25024	PRX27024	PRX30024	PRX35024	PRX40024
26	PRX15026	PRX17026	PRX20026	PRX25026	PRX27026	PRX30026	PRX35026	PRX40026
30	PRX15030	PRX17030	PRX20030	PRX25030	PRX27030	PRX30030	PRX35030	PRX40030
31	PRX15031	PRX17031	PRX20031	PRX25031	PRX27031	PRX30031	PRX35031	PRX40031
32	PRX15032	PRX17032	PRX20032	PRX25032	PRX27032	PRX30032	PRX35032	PRX40032
34	PRX15034	PRX17034	PRX20034	PRX25034	PRX27034	PRX30034	PRX35034	PRX40034
38	PRX15038	PRX17038	PRX20038	PRX25038	PRX27038	PRX30038	PRX35038	PRX40038