



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610 Fax: 86-532-5568 9611 E-mail: med@allpro.com.cn



www.allprocorporation.com



ABSORBENT GAUZE BANDAGE

Description:

1. Material: 100% cotton
 2. Width: 2.5cm, 5cm, 7.5cm, 10cm, 15cm, 20cm
 3. Length: 4.5m, 5m, 10m
 4. Density: 40's*40's; 19*10mesh, 19*15mesh, 26*18mesh, etc
Other dimensions, width, length, and package can be produced as per requirement
 5. Description: A gauze bandage is a thin, woven fabric material that is placed over a wound to keep it clean while allowing air to penetrate and promote healing. It can be used to secure a dressing in place, or it can be used directly on a wound. These bandages are the most common type and are available in many sizes. They are good for hard to reach areas that may be difficult to dress.
-