



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610 Fax: 86-532-5568 9611 E-mail: med@allpro.com.cn



www.allprocorporation.com



Description:

3*6cm pad with 5*5cm pouch, 70% isopropyl alcohol, individually sealed swabs, 100pcs/box, 200 pcs/box, 10000pcs/carton
Alcohol swab, with 70% IPA.
Original solution from Exxon Mobil 99.8% solution.
60x30mm, single pack.

Standard packing:

English neutral packing, 100pc/box or 200pc/box,
100box/carton;
or as client requirement.

Specs:

Pouch size: 50*50±2mm
Pouch Material: Aluminum foil paper, 73g/m²
Non-woven size: 60*30±2mm,
Solution: 70% isopropyl alcohol
Solution Weight: 0.4g- 0.5g
The material of the box: cardboard with white face
and mottled back; 300g/m²;