



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610 Fax: 86-532-5568 9611 E-mail: med@allpro.com.cn



www.allprocorporation.com



## ABSORBENT GAUZE ROLL

### Description:

1. Material: 100% cotton
2. Accessories: X-ray detectable or without
3. Width: 90cm-4ply, 120cm-4ply
4. Length: 50m, 100m, 200m, 500m
5. Density: 40's\*40's; 19\*10mesh, 19\*15mesh, 26\*18mesh, etc  
Other dimensions, width, length, and package can be produced as per requirement
6. Description: The gauze is degreased and bleached by advanced way to ensure superior purity and absorbency. The quality meet the standard of English Medical Dictionary. The product has no fluorescence .It is used widely in medical circle and other area.