



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Ryle's Stomach Tube



Standard Type

STOMACHE TUBE RYLE'S STOMACH TUBE /STANDARD TYPE

Description:

- X- ray line/ radio-opaque line is available for accurate placement
- Tube end is sealed with stainless steel balls, to assist the catheter indwelling in a certain place

Item No.	Item No. (DEHP FREE)	Size (Fr)	Color Code
APF-1	F 250 1200	12	White
APF-2	F 250 1400	14	Green
APF-3	F 250 1600	16	Orange
APF-4	F 250 1800	18	Red
APF-5	F 250 2000	20	Yellow
APF-6	F 250 2200	22	Blue
APF-7	F 250 2400	24	Violet