



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



THORACIC DRAINAGE TUBE

Description:

- Non-toxic, medical grade PVC material
- X-ray opaque line.
- Type: round or fish-mouth tip
- Length: 500mm
- Size: 20/22/24/26/28/30/32/34/36