



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH SMALL WINGS

Description:

F. E. P material

Thin-wall construction permits high flow rates

Radiopaque, x-ray detectable stripes

Clear and transparent

Smooth surfaces

Back-cut special designed needle, siliconized

Rotated wing for fixation and as grip plate

Hydrophobic blood stopper

Color coded

Lock fitting

Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G

Package: 30pcs/box, 600pcs/ctn, 77*53*29cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA Y TYPE

Description:

- Thin-wall construction permits high flow rates
- Radiopaque, x-ray detectable stripes
- Clear and transparent
- Smooth surfaces
- Back-cut special designed needle
- Hydrophobic blood stopper
- Color coded
- Lock fitting
- Plastic packing and blister packing are provided
- Reduce insertion attempts
- Increase protection against blood exposure and needle stick injuries
- Simplify IV therapy process
- Enhance patient safety and satisfaction

Size: 14G/16G/17G/18G/20G/22G/24G

Package: 50pcs/box, 1000pcs/ctn, 50.5*42.5*28cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH INJECTION PORT

Description:

F.E.P material
Thin-wall construction permits high flow rates
Radiopaque, x-ray detectable stripes
Clear and transparent
Smooth surfaces
Back-cut special designed needle, siliconized
Flexible fixation wings
Injection valve with cap
Hydrophobic blood stopper
Color coded
Lock fitting
Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G

Package: 50pcs/box, 1000pcs/ctn, 48*42*32cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
14G	Orange	2.10	45	240
16G	Grey	1.74	45	180
17G	White	1.40	45	125
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA PEN LIKE

Description:

F.E.P material

Thin-wall construction permits high flow rates

Radiopaque, x-ray detectable stripes

Clear and transparent

Smooth surfaces

Back-cut special designed needle, siliconized

Hydrophobic blood stopper

Color coded

Lock fitting

Sterilization: Ethylene Oxide Gas.

Size: 14G/16G/17G/18G/20G/22G/24G

Package: 50pcs/box, 1000pcs/ctn, 50*34*33cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V. CANNULA SAFETY TYPE

Resoted in the body as combo with infusion sets for transfusion.

This product can be kept in the vein of human for a long time.

it brings less vascular injury. Consequently it largely reduces suffering of patients and reduces intensity of health care work.

The needle retract into the tube after single use, and prevent reusing and protect medical care personnel

Size: 16G/18G/20G/22G/24G


Package:50pcs/box,1000pcs/ctn



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V. CANNULA SAFETY CLOSED LIGHT PROOF

Resoted in the body as combo with infusion sets for transfusion.

This product can be kept in the vein of human for a long time.

It brings less vascular injury. Consequently it largely reduces suffering of patients and reduces intensity of health care work.

The needle retract into the tube after single use, and prevent reusing and protect medical care personnel
Closed light proof to infuse.



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D. mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA PEN-LIKE SAFETY MODEL

Thin-wall construction permits high flow rates

Radiopaque, 3x-ray detectable stripes

Clear and transparent

Smooth Surfaces

Back-cut special designed needle

Color coded

Lock fitting

Plastic packing and Blister packing are provided

Locking function avoid secondary use, protect the harm of

Size: 16G/18G/20G/22G/24G

Package: 50pcs/box, 1000pcs/ctn, 48.5*30.5*37.5cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH WINGES SAFETY TYPE

- Thin-wall construction permits high flow rates
- Radiopaque,3x-ray detectable stripes
- Clear and transparent
- Smooth Surfaces
- Back-cut special designed needle
- Rotated wing for fixation and as grip plater
- Color coded
- Lock fitting
- Plastic packing and Blister packing are provided
- Locking function avoid secondary use,protect the harm of

Size: 16G/18G/20G/22G/24G

Package:50pcs/box, 1000pcs/ctn, 60*28.5*52cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Gauge	Colour code	Catheter O.D.mm	Catheter Length mm	Water Flow Rate(ml/min)
16G	Grey	1.74	45	180
18G	Green	1.20	45	80
20G	Pink	1.11	32	54
22G	Blue	0.80	25	31
24G	Yellow	0.60	19	19
26G	Dark Blue	0.60	19	13

I.V.CANNULA WITH INJECTION VALVE SAFETY TYPE

Thin-wall construction permits high flow rates

Radiopaque, 3x-ray detectable stripes

Clear and transparent

Smooth Surfaces

Back-cut special designed needle

Flexible fixation wings

Injection valve with cap

Color coded

Lock fitting

Plastic packing and Blister packing are provided

Locking function avoid secondary use, protect the harm of

Size: 16G/18G/20G/22G/24G

Package: 50pcs/box, 1000pcs/ctn, 41*32*47cm